



# REDUCE LAP-HF II

Corvia® Atrial Shunt in heart failure with preserved or mildly reduced ejection fraction

## 2-YEAR RESPONDER GROUP HIGHLIGHTS<sup>1</sup>

### 50%

Reduction in HF event rate compared to sham control

### 46%

Greater improvement in KCCQ-OSS\* compared to sham control

### 98%

Patency



### CONCLUSION

The Corvia Atrial Shunt improved quality of life and reduced heart failure (HF) symptoms and events through 2 years without a safety signal in appropriately selected HF patients with an ejection fraction (EF)  $\geq 40\%$ .



### BACKGROUND

REDUCE LAP-HF II is the world's first phase III randomized trial for atrial shunting in HF. It defined the optimal patient population for atrial shunt therapy (the "responder" population) and demonstrated the Corvia Atrial Shunt was safe, significantly reduced HF events, and improved quality of life in the responder population compared to sham control at one year.<sup>2</sup>




### RELEVANCE

This two-year data, the longest follow-up from a randomized controlled trial on atrial shunt therapy, provides valuable insight into the Corvia Atrial Shunt's long-term potential to meet the significant unmet need in heart failure patients with EF  $\geq 40\%$ .

### STUDY POPULATION

Randomized, double-blind, placebo-controlled  
**621 patients**



Atrial Shunt




Sham Control

**PATIENTS INCLUDED:**  $\geq 40$  years with symptomatic HF and EF  $\geq 40\%$  with elevated exercise PCWP ( $\geq 25$  mm Hg) and left-to-right gradient ( $\geq 5$  mm Hg)

#### RESPONDER COHORT Identified (n=313)

Patients without significant pulmonary vascular disease (PVR  $< 1.74$  WU) and no pacemaker





Overall Population

← Follow-Up Blinding Maintained Through 2 Years →

\*Kansas City Cardiomyopathy Questionnaire Overall Summary Score


# Responder Population Outcomes Through 2 Years<sup>1,2</sup> (n=313)

## HEART FAILURE EVENTS



Shunt therapy led to a 50% reduction in the rate of HF events


## QUALITY OF LIFE



Shunt patients sustained improved health status with 46% greater improvement (+4.6 points) in KCCQ

## Learn More

Prof. Finn Gustafsson and Prof. Maja Cikes share their thoughts on the 2-year data from REDUCE LAP-HF II



## KEY SAFETY OUTCOMES

Events in Responders through 2 years	Corvia Atrial Shunt (N=161)	Sham control (N=152)	p-value
Cardiovascular mortality	3/160 (1.9%)	2/148 (1.4%)	0.72
Non-fatal ischemic stroke	3/160 (1.9%)	0/148 (0.0%)	--
New or worsening kidney dysfunction	14/160 (8.8%)	22/148 (14.9%)	0.10
Major adverse cardiac events	6/160 (3.8%)	4/148 (2.7%)	0.61
Thrombo-embolic complications	1/160 (0.6%)	1/148 (0.7%)	0.96
≥30% Decrease in TAPSE	5/160 (3.1%)	4/148 (2.7%)	0.83



Corvia Medical, Inc.  
One Highwood Drive, Suite 300  
Tewksbury, MA 01876 USA

IASD and Corvia are registered trademarks of Corvia Medical, Inc.  
©2024 Corvia Medical, Inc. All rights reserved.

PS00698, MKT1020 (US/CA/AU) Rev01 2024-08

- Gustafsson, F, Petrie, M, Komtebedde, J. et al. 2-Year Outcomes of an Atrial Shunt Device in HFpEF/HFmrEF: Results From REDUCE LAP-HF II. *J Am Coll Cardiol HF*. 1 Jun 2024. <https://doi.org/10.1016/j.jchf.2024.04.011>
- Borlaug, BA, Blair, J, Bergmann, MW et al. Latent Pulmonary Vascular Disease May Alter the Response to Therapeutic Atrial Shunt Device in Heart Failure. *Circulation*. 2022;10.1161/CIRCULATIONAHA.122.059486.
- Incidence rate ratio (IRR) at 1- and 2-years = 0.49 [95% CI: 0.25-0.95], p=0.034
- Statistical analyses conducted by Baim Institute for Clinical Research. Data on file.

CAUTION: Investigational Device. Limited by United States law to investigational use. To be used by qualified investigators only. For use in a pre-market clinical investigation only.